



Genlab Drying/Warming Cabinets Power Rating

Drying/Warming Cabinets (Non Insulated)			
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>	
DC/125	125	500	
DC/250	250	750	
DC/500	500	1750	
DC/1000	1000	2500	

Drying/Warming Cabinets (Insulated)			
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>	
IWC/100	100	500	
IWC/200	200	750	
IWC/425	425	1750	
IWC/885	885	2500	

***Important Note**

Due to the enhanced thermal efficiency of the insulated drying cabinet will consume significantly less energy during the working life of the unit.

