



Genlab General Purpose Incubators Power Ratings

Vertical style		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
Mini/6	6	150
Mini/18	18	350
Mini/30	30	500
Mini/40	40	500
Mini/50	50	500
Mini/75	75	500



Horizontal style		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
INC/50	50	500
INC/75	75	500
INC/100	100	500
INC/125	125	750
INC/150	150	750
INC/200/F	200	1000



***Important Note**

All Genlab ovens and incubators are designed and built with the environment in mind. Once operating temperature has been reached, the oven/incubator only consumes power for approximately 20% of the running time, due to their thermal efficiency.

