



## Genlab Large Capacity Multi – Purpose Incubators Power Rating

<b>Large Capacity Multi – Purpose Incubators (Horizontal Style)</b>		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
LCI/9H/DIG	250	1100
LCI/12H/DIG	350	1500
LCI/15H/DIG	425	2000
LCI/18H/DIG	500	2000
LCI/22H/DIG	650	2000
LCI/27H/DIG	750	2000
LCI/30H/DIG	850	2000
LCI/35H/DIG	1000	2000
LCI/42H/DIG	1250	2000



<b>Large Capacity Multi – Purpose Incubators (Vertical Style)</b>		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
LCI/9V/DIG	250	1100
LCI/12V/DIG	350	1500
LCI/15V/DIG	425	2000
LCI/18V/DIG	500	2000
LCI/22V/DIG	650	2000
LCI/27V/DIG	750	2000
LCI/30V/DIG	850	2000
LCI/35V/DIG	1000	2000



### **\*Important Note**

All Genlab ovens and incubators are designed and built with the environment in mind. Once operating temperature has been reached, the oven/incubator only consumes power for approximately 20% of the running time, due to their thermal efficiency.

