



Genlab Large Capacity Multi – Purpose Ovens Power Rating

Large Capacity Multi – Purpose Ovens (Horizontal Style)		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
LCO/9H/DIG	250	2000
LCO/12H/DIG	350	2500
LCO/15H/DIG	425	3000
LCO/18H/DIG	500	3000
LCO/22H/DIG	650	4500
LCO/27H/DIG	750	6000
LCO/30H/DIG	850	6000
LCO/35H/DIG	1000	6000
LCO/42H/DIG	1250	6000



Large Capacity Multi – Purpose Ovens (Vertical Style)		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
LCO /9V/DIG	250	2000
LCO /12V/DIG	350	2500
LCO /15V/DIG	425	3000
LCO /18V/DIG	500	3000
LCO/22V/DIG	650	4500
LCO/27V/DIG	750	6000
LCO/30V/DIG	850	6000
LCO/35V/DIG	1000	6000



***Important Note**

All Genlab ovens and incubators are designed and built with the environment in mind. Once operating temperature has been reached, the oven/incubator only consumes power for approximately 20% of the running time, due to their thermal efficiency.

