



## Genlab Cooled Incubators/Water Baths Power Ratings

<b>Cooled Incubators</b>			
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>	
Ci/50	50	300	
CI/100	100	300	
CI/200	200	400	
Ci/280	280	550	

<b>Water Baths</b>			
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)</i>	
WBH6	6	750	
WBH9	9	750	
WBH14	14	750	
WBH22	22	1500	

### **\*Important Note**

All Genlab ovens and incubators are designed and built with the environment in mind. Once operating temperature has been reached, the oven/incubator only consumes power for approximately 20% of the running time, due to their thermal efficiency.

