

# Genlab

## Laboratory Division

### COOLED INCUBATORS



50 litre Cooled Incubator  
CI/50

100 litre Cooled Incubator  
CI/100

The GENLAB Cooled Incubators offer a range of highly accurate, efficient and economical units that are ideally suited for a wide range of research, testing and general laboratory applications.

#### Features

- |  |   |
|--|---|
| ■ 4 sizes: 50, 100, 200 & 280 litre  | ■ Self evaporating condensate drip-tray |
| ■ Temperature range +5°C to 40°C   | ■ RF and TV suppressed                  |
| ■ Fluctuation +/- 0.1°C @ 20°C   | ■ Sealed refrigeration system           |
| ■ PID Microprocessor digital control with dual display of set point and actual temperature | ■ Automatic defrost system              |
| ■ Auto-switch for refrigeration  | ■ Locking door                          |
| ■ White stove enamel exterior  | ■ Easy clean plastic moulded chamber    |
| ■ Fan assisted circulation   | ■ Safety overheat thermostat            |
|  | ■ One year warranty                     |
|  | ■ C E compliant                         |

#### Accessories

- |                                  |                              |
|----------------------------------|------------------------------|
| ■ Timers                         | ■ Internal 13 amp socket     |
| ■ Fully programmable controllers | ■ Audible and visible alarms |
| ■ Communication ports            | ■ Viewing window             |
| ■ Interior lights                | ■ Chart recorders            |
| ■ Access ports                   |                              |



Certificate No. 5531