

Genlab

Laboratory Division

COOLED INCUBATORS



50 litre Cooled Incubator
CI/50



100 litre Cooled Incubator
CI/100

The GENLAB Cooled Incubators offer a range of highly accurate, efficient and economical units that are ideally suited for a wide range of research, testing and general laboratory applications.

Features

- | | |
|--|---|
| ■ 4 sizes: 50, 100, 200 & 280 litre | ■ Self evaporating condensate drip-tray |
| ■ Temperature range +5°C to 40°C | ■ RF and TV suppressed |
| ■ Fluctuation +/- 0.1°C @ 20°C | ■ Sealed refrigeration system |
| ■ PID Microprocessor digital control with dual display of set point and actual temperature | ■ Automatic defrost system |
| ■ Auto-switch for refrigeration | ■ Locking door |
| ■ White stove enamel exterior | ■ Easy clean plastic moulded chamber |
| ■ Fan assisted circulation | ■ Safety overheat thermostat |
| | ■ One year warranty |
| | ■ C E compliant |

Accessories

- | | |
|----------------------------------|------------------------------|
| ■ Timers | ■ Internal 13 amp socket |
| ■ Fully programmable controllers | ■ Audible and visible alarms |
| ■ Communication ports | ■ Viewing window |
| ■ Interior lights | ■ Chart recorders |
| ■ Access ports | |



Certificate No. 5531