



Genlab General Purpose Ovens Power Ratings

Vertical style		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
Mino/6	6	350
Mino/18	18	500
Mino/30	30	750
Mino/40	40	750
Mino/50	50	750
Mino/75	75	1000



Horizontal style		
<i>Genlab Ref</i>	Capacity (litres)	<i>Power Rating (W)*</i>
OV/50	50	750
OV/75	75	1000
OV/100	100	1000
OV/125	125	1250
OV/150	150	1500
OV/200/F	200	2000



***Important Note**

All Genlab ovens and incubators are designed and built with the environment in mind. Once operating temperature has been reached, the oven/incubator only consumes power for approximately 20% of the running time, due to their thermal efficiency.

