




Data Logging Software

# G - Viewer




The Genlab G-Viewer allows you to display oven data into both graphical and numeric presentation form. Designed bespoke for the Genlab range of Incubators, Ovens and Drying Cabinets, it is intuitive to use by downloading data directly from Genlab product ranges via USB.

As well as process data, the set point data, energy use and audit data are presented. With user selectable time lines, the data can show minimum, maximum, average and energy used per process.

## Features

- Temperature, Audit and Energy Used Displayed
- Intuitive to use
- Licence Free
- Exportable to Excel
- Batch Data Saveable



Doc: GViewer Rev 1

Tanhouse Lane, Widnes, Cheshire WA8 0SR  
Tel 0151 424 5001 enquiries@genlab.co.uk www.genlab.co.uk

