




Front View



Side View



Rear



Top View

- Notes
- 1 - Dimensions shown in mm unless otherwise stated
 - 2 - Supply - 230V
 - 3 - Load - 0.25A
 - 4 - Airflow - 30.5 litres / sec

This Drawing is the property of GENLAB Limited and cannot be copied or used without their written consent					
	Title:	Genlab Extractor	Drawn by:	Checked by:	Approved by:
	Client:	Stock Item	Rick J		
	Drawing No:	GL-EXT	Rev:	Date:	Sheet No.
	Contract:	Rick	1	27/07/21	1 of 1