


Materials Testing Division

Riffle Boxes



These sample dividers are constructed of heavy gauge sheet metal finished in an easy clean powder coated paint and are supplied in a range to meet the requirements of BS 1377, 1924, 812 and ASTM C136.

Each riffle is supplied complete with divider and three containers.

Genlab ref. Number	Slot width (mm)	No. of slots	Genlab ref. Number	Slot width (mm)	No. of slots
RB607	7	12	RB630	30	10
RB610	10	12	RB638	38	8
RB613	13	12	RB640	40	8
RB615	15	12	RB645	45	8
RB620	20	10	RB650	50	8
RB625	25	10	RB665	65	8